


Indel-based Realignment

Improving the original alignments of the
reads based on multiple sequence
(re-)alignment

We are here in the Best Practices workflow

REALIGNMENT

Calling Variants with the GATK 2




* New tools or functionalities not available in GATK-Lite

PURPOSE


Why realign around indels?

- InDels in reads (especially near the ends) can trick the mappers into mis-aligning with mismatches
 - These artifactual mismatches can harm base quality recalibration and variant detection (unless a sophisticated caller like the Haplotype Caller is used)
- ☒ **Realignment around indels helps improve the accuracy of several of the downstream processing steps.**

An example of a strand-discordant locus



Local realignment uncovers the hidden indel in these reads and eliminates all the potential FP SNPs




PRINCIPLES

Three types of realignment targets

- Known sites (e.g. dbSNP, 1000 Genomes)
- Indels seen in original alignments (in CIGARs)
- Sites where evidence suggests a hidden indel

Local realignment identifies most parsimonious alignment along all reads at a problematic locus

1. Find the best alternate consensus sequence that, together with the reference, best fits the reads in a pile (maximum of 1 indel)




2. The score for an alternate consensus is the total sum of the quality scores of mismatching bases
3. If the score of the best alternate consensus is sufficiently better than the original alignments (using a LOD score), then we accept the proposed realignment of the reads

PROTOCOL

Indel Realignment steps/tools


- Identify what regions need to be realigned

→ **RealignerTargetCreator**




- Perform the actual realignment


→ **IndelRealigner**



Indel Realignment workflow



Indel Realignment workflow




RealignerTargetCreator

- Pre-processing step to find intervals that may need realignment


```
java -jar GenomeAnalysisTK.jar -T RealignerTargetCreator \  
  -R human.fasta \  
  -I original.bam \  
  -known indels.vcf \  
  -o realigner.intervals
```

- Input BAM file not necessary if processing only at known indels
- Using a list of known indels will both speed up processing and improve accuracy, but is not required

Indel Realignment workflow



Indel Realignment workflow




IndelRealigner

- Attempts realignment at RTC target intervals

```
java -jar GenomeAnalysisTK.jar -T IndelRealigner \  
    -R human.fasta \  
    -I original.bam \  
    -known indels.vcf \  
    -targetIntervals realigner.intervals \  
    -o realigned.bam
```

- Must use same input file(s) used in RealignerTargetCreator step
- Processing options
 - Only at known indels: much faster, accurate for ~90-95% of indels
 - At indels seen in the original BAM alignments: the recommended mode
 - Using full Smith-Waterman realignment: most accurate, but heavy computational cost and not really necessary with the new techs

Indel Realignment workflow




RESULTS



This is what a realigned BAM looks like

NA12878, chr1:1,510,530-1,510,589

Before



After



Did the realignment work properly?


- Indel Realigner changes the CIGAR string of realigned reads but maintains the original CIGAR (with OC tag)
 - So it's very easy to check that realignment was performed and/or how many reads were adjusted
- BUT no formal measure to assess the accuracy or completeness of the realignment process

Is realignment still necessary with latest software?

- Latest tools being implemented for discovering mutations all include some sort of assembly step (for which upstream realignment is not really helpful).
- BUT big improvement for Base Quality Score Recalibration when run on realigned BAM files (artifactual SNPs are replaced with real indels).
- Also still useful for legacy tools, e.g. full realignment should be performed if using the GATK's Unified Genotyper.

We were here in the Best Practices workflow

NEXT STEP: BASE RECALIBRATION



Further reading

<http://www.broadinstitute.org/gatk/guide/topic?name=intro>

<http://www.broadinstitute.org/gatk/guide/topic?name=best-practices>

<http://www.broadinstitute.org/gatk/guide/article?id=38>

[http://www.broadinstitute.org/gatk/gatkdocs/
org_broadinstitute_sting_gatk_walkers_indels_IndelRealigner.html](http://www.broadinstitute.org/gatk/gatkdocs/org_broadinstitute_sting_gatk_walkers_indels_IndelRealigner.html)

[http://www.broadinstitute.org/gatk/gatkdocs/
org_broadinstitute_sting_gatk_walkers_indels_RealignerTargetCreator.html](http://www.broadinstitute.org/gatk/gatkdocs/org_broadinstitute_sting_gatk_walkers_indels_RealignerTargetCreator.html)